


BURN S

Burn stop dressing



The natural plant extract gel absorbs the heat and cools the injured zone, no matter how severe the damage is, it will stop the burning process of the wound caused by the heat

Gel basically made of water with different physical characteristics: the humidity prevents the gauze being attached to the wound making the removal of the gauze extremely easy and painless any time you need to inspect it or cure it

Special soft and elastic support, saturated with gel that avoids listelling

By applying the gel you will obtain a very efficient cooling effect on the burnt surface and it guarantees:

- reduction of the deepness of the burn
- reduction of the aedema
- reduction of pain

The sterile package offers the advantage of an efficient protection against infections and as a medication kit.

The Burn S dressing is a new burn stop idrogel made of natural plant extracts. A sterile single packed dressing, easy to keep and always ready to use.

TECHNICAL DATA

Burn S (BA) 20 x 20 cm [7.87 x 7.87 in]

BS04021D

