

BURN C

Fire blanket



This blanket can be used as protection for rescue workers entering or leaving a fire area. It can also be used to extinguish flames on a burn victim and, if necessary, be used as a large dressing, providing instant pain relief and beginning the healing process while victims are being transported to emergency care units.

Specific Characteristics

- Sterile trauma blanket
- Shoulder strap make it easy to carry
- Suitable for: extinguish flame, large dressing, entry or exit from fire
- No loose fibers
- Relieves pain
- Cools and provides the necessary water supply to the skin
- Non-adherent and non-toxic

Technical data

Dimensions	245x160 cm
Idrogel composition	96% water, Melaleuca Alternifolia, emulsifiers
PH Gel	5,5-7
Sterilized by	Gamma ray
Blanket Material	Wool

Class IIb MD according to Dir. 93/42/CEE

BURN C 160 X 254 CM - FIRE BLANKET
BS04040



CND Classification	M040499
Registration number	67954

Distributed by Spencer Italia S.r.l.

Rev.1 (26/10/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it