

OXYMASK 150

High concentration oxygen masks



Disposable face mask for oxygen therapy.

Specific Characteristics

- Adjustable elastic
- Semi transparent shell to improve patient monitoring
- Anti-crush tube
- Flexible aluminum nose clip to increase stability
- Comes with reservoir

Latex Free

Technical data

Material	PVC
Tube length	2 m

Class IIa Med.Device according to Dir. 93/42/CEE

ADULT OXYGEN MASK HIGH
CONCENTRATION WITH RESERVO
OX20104A



CND Classification	R03010201
Registration number	104794

Rev.1 (14/02/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it