

B-BAG

Double bag resuscitator



Double bag self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The double bag transmits the pressure variations to the operator more effectively
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Steam sterilizable

Technical data

Materials	Rubber, Silicone, PC, PVC
Bag volume	1600 ml
Reservoir volume	2600 ml
Weight of the kit	490 g

Standard equipment

BG07002	Res.bag w/valve for B-Bag adult
EV50011	EVX 11 – Non reb. Overpr. valve

Med. Dev. CE marked according to Dir. 93/42/CEE

B-BAG - ADULT DOUBLE BAG
 W/SOVR. RESERV.
BG07020B



Registration number 58903

Optionals

OX20131A	OXX30 anti crushing oxygen tube 10pcs
----------	---------------------------------------

Rev.0 (07/04/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it