

B-BAG

Double bag resuscitator



Double bag self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The double bag transmits the pressure variations to the operator more effectively
Possibility to connect the bag to sources of oxygen supply
Low storage space required

Technical data

Materials	Rubber, Silicone, PC, PE
Bag volume	280 ml
Overpressure valve	40 cmH ₂ O
Weight of the kit	200 g

Standard equipment

EV50007 | EVX 07 – Non reb. Overpr. valve

Med. Dev. CE marked according to Dir. 93/42/CEE

B-BAG NEONATAL MANUAL
RESUSCITATOR W/OVERPR.VALV
BG07021B



Registration number | 58904

Optionals

OX20131A | OXX30 anti crushing oxygen tube 10pcs

Rev.1 (07/04/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it