

## B-Life

### Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands
Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications
Connection standard ISO 15/22 mm
Possibility to connect the bag to sources of oxygen supply
Comes with non-rebreathing valve
Low storage space required
Steam sterilizable

#### Technical data

Materials	Silicone, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Stroke volume range	300 ml with 1 hand 350 ml with 2 hands
Maximum oxygen concentration <sup>(with reservoir)</sup>	98%
Max. re expansion rhythm	100 bpm

#### Optionals

BG03091A	Reservoir
RM20403B	B-LIFE – Silicone face Mask sz.3
OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
EV50007C	EVX 07 – Overp. and non-rebreathing valve

MD CE marked according to Dir. 93/42/CEE

B- LIFE SILICON PAEDIATRIC  
RESUSCITATION BAG  
BG03004A



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)