

B-Life

Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use INFANT marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Low storage space required
 Latex Free
 Steam sterilizable

Standard equipment

EV50007C EVX 07 - Overpressure and non-rebreathing valve

Optionals

OX20131A OXX30 anti crush oxygen tubes 10pcs
 BG02001A Laxan transparent bag for manual resuscitator

Technical data

Materials	Silicone, PC
Bag Volume	280 ml
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration	95%
Max. re expansion rhythm	95 bpm

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE INFANT
 W/OVERP. VALVE R-VALVE
 BG03025A



Registration number 34003

Rev.1 (07/06/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it