

## B-Life

### Pediatric manual resuscitation unit



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

- The material of construction and the surface finish increase grip to the operator's hands
- Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications
- Connection standard ISO 15/22 mm
- Possibility to connect the bag to sources of oxygen supply
- Autoclavable (except reservoir)
- Low storage space required

#### Standard equipment

BG03091A	Reservoir
EV50007C	EVX 07 - Overpressure and non-rebreathing valve
RM20403B	B-LIFE – Silicone face Mask sz.3

#### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

#### Technical data

Materials	Silicone, Medical grade PVC, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Stroke volume range	300 ml with 1 hand 350 ml with 2 hands
Overpressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration <small>(with reservoir)</small>	98%
Max. re expansion rhythm	100 bpm
Overall weight	290 ±10 g

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE  
PAED.W/OVERP.RESERV. R-VALVE  
AND MASK SIZE 3



**BG03033B**

Registration number 21706

Rev.3 (29/09/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)