

## B-Life

### Neonatal manual resuscitation unit



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The material of construction and the surface finish increase grip to the operator's hands .  
 Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications  
 Connection standard ISO 15/22 mm  
 Possibility to connect the bag to sources of oxygen supply  
 Low storage space required

#### Technical data

Materials	Silicone, medical grade PVC, rubber , PC
Bag volume	280 ml
Reservoir volume	600 ml
Over pressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration <sup>(with reservoir)</sup>	95%
Max. re expansion rhythm	95 bpm
Overall weight	310±30 g

#### Standard equipment

BG03092A	Reservoir
EV50007C	EVX 07 – Overpressure and non-rebreathing valve
RM20504A	ECO B-LIFE – Rubber face mask sz.0
RM20503A	ECO B-LIFE – Rubber face mask sz.1
RM20502A	ECO B-LIFE – Rubber face mask sz.2
OX20130A	Anti-crushing oxygen tube

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SIL.RESUSCITATOR INFANT  
 W/OVERPR. RES. MASK OXYGEN  
 TUBE R-VALVE



**BG03044A**

#### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.0 (20/02/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)