

PRO B-Life

PVC Resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Latex Free

Technical data

Materials	Medical grade PVC, PC
Bag Volume	1600 ml
Stroke volume	900 ml
Overpressure valve	60 cmH ₂ O
Max. re expansion rhythm	45 bpm

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator
BG03091A	Reservoir bag
RM10894A	Cuffed mask Size 5

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC ADULT
 RESUSCITATOR WITH R-VALVE
BG02006A



Registration number 28803

Rev.0 (04/09/2013)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it