

# Lienz Krone

Reinforced aluminium case for resuscitation equipment complete adult



## Krone Series

They are characterized by:

- supports for improved stability
- three hooks for mounting on vehicles or for attachment of shoulder belts
- reinforced angles in stainless steel
- robust and comfortable handles
- zip protection
- two key-locks
- mobile thermo-formed service tray
- phial holder in shockproof rubber
- adjustable shoulder belt with shoulder protector
- sticker kit for personalization

## Lienz Krone

- Transparent shockproof panel
- Mobile dividing shockproof panels
- Straps for fixing the oxygen tank

Contents:

Cylinder 2 L with reducer and flow selector  
 Oxygen therapy mask Adult  
 Pro B-Life resuscitator PVC adult  
 Reanimation mask Air Cuffed Mask size 3  
 Guedel size 1  
 Guedel size 2  
 Guedel size 3  
 Collin tongue holder  
 Plastic tongue stick  
 Traditional laryngoscope handle  
 Curved blade size 1  
 Curved blade size 2  
 Curved blade size 3  
 Magill pince adult  
 Basic Bp phonendoscope  
 Basic Bp sphygmomanometer  
 Latex tubular tourniquet  
 Isothermic sheet adult  
 Universal scissors  
 Pen light  
 Hydrogel burn dressing 10 x 10 cm  
 Res-Q.Vac suction unit  
 Pressure flow selector

Dimensions	520 x 225 x h370 mm
Weight complete	13,5 kg
Capacity	42,3 L
Material	Aluminium

Lienz Krone CB06018A



Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

E-mail: [info@spencer.it](mailto:info@spencer.it) Web: <http://www.spencer.it>

This document has been automatically generated by [www.spencer.it](http://www.spencer.it) on 2010/1/13 12:12:19, and is referred to the last production lot. The informations furnished have been last updated on 2009/12/01. To verify the correctness of this document, please contact our Quality Dept. at [quality@spencer.it](mailto:quality@spencer.it)