

## B-BAK PIN MAX

Spine board complete of pins



A special B-bak, very much appreciated by rescue teams, it puts together its famous technical characteristics and a built in system that allows to fix the restraints to the board.

The side handles are translucent incorporated elements, they allow you easy buckling-unbuckling operations for a fast and safe anchorage of the patient's restraints with shackle.

### OPTIONAL ACCESSORIES

Personalising System.

ST02106 MARK-UP

Restraint belts for pin anchorage

ST02015 PIN STRAPS

### TECHNICAL DATA

Material	HDPE
Dimensions	405 x 45 x h1840 mm [15.94 x 1.96 x h72.44 in]
Weight:	6 kg [13.22 lbs]
Load capacity B-bak	
Pin MAX	454 kg [1000 lbs]
Anchor Pins	4
Number of handles	14
X ray transparency	Yes
Functioning temperature	from -10 to +50 °C [from -14 to +122 °F]

B-bak Pin MAX yellow

ST02067B

