

## B-BAK PIN MAX MILITARY

Maximum loading weight spine board with pins



The special B-bak, a favourite of the rescue squads which is a combination of the renowned technical characteristics for which the product is renowned and an integrated fixing system for the restraint straps.

The inbuilt side supports are made with a x ray compatible material and permit quick and safe attaching and un attaching of the patient restraint belts thanks to the carabineer end hook.

### OPTIONAL ACCESSORIES

Personalising System.

Restraint belts for pin anchorage

ST02106 MARK-UP

ST02015 PIN STRAPS

### TECHNICAL DATA

Material	Polyethylene
Dimensions	405 x 45 x h1840 mm [15.94 x 1.96 x h72.44 in]
Weight:	6 kg [13.22 lbs]
Load capacity B-bak	
Pin MAX	454 kg [1000 lbs]
Anchor Pins	4
Number of handles	14
X ray transparency	Yes
Functioning temperature	from -10 to +50 °C [from 14 to +122 °F]

B-BAK PIN MAX MILITARY ST02068B

