

# EVA

## Resuscitation kit



Compact semi-rigid professional backpack equipped with devices most commonly used in resuscitation procedures.

Included Equipment	
BG01009	ADULT RESUSCITATOR
BG01010	PEDIATRIC RESUSCITATOR
EV50010	N/REBREATHING VALVE
EV50009	ADULT OVERPRESSURE VALVE
RM10883	CUFFED MASK SZ 3
RM10885	CUFFED MASK SZ 5

Specific Features	
Semi-rigid case to provide more protection to the content	
Backpack with 2 adjustable shoulder straps	
Internal separator with elastic bands	
Mash pocket on the internal side of the front half of the backpack	
Handle on the upper side for an easy transport	
Technical data	
Dimensions	380x300x200 ± 30 mm
Materials	EVA, Nylon
Weight of empty backpack	0,73 ± 0,1 kg

Medical devices CE marked according to dir. 93/42/CE	
EVA BLUE RESUSCITATION KIT CO B-LIFE - EMERGENCY BACKPACK	
<b>CB00216A</b>	

Rev.0 (18/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)