

## QMX 120

Pump for vacuum systems



Pump intended for use with with Spencer vacuum mattresses or splints.

### Specific features

The aluminum chamber makes it a rugged and lightweight device

The shape of the base allows to hold the pump firmly during use

### Technical data

Dimensions	500x140x80 mm
Tube length	100 ± 5 cm
Materials	Al, PVC, Steel
Weight	1,14 ± 0,1 kg

ALUMINUM VACUUM PUMP  
QMX120  
QM22120A

Rev.1 (05/12/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)