

**QMX 125**

Double-act vacuum pump



Double acting pump that allows to depressurize two devices at the same time. Aluminium body, zinc laminated anti-slip pedal, plunger mechanically fixed to the knob.

Characteristics that permit the necessary strength to generate void in vacuum mattresses.

The connecting tube with a wide gauge, attached to an unique one-way valve, guarantees good flow and speeds up rescue operations.

The building technique, the construction materials and the attention to details guarantee a long-life span.

**TECHNICAL DATA**

Height	500 mm
Weight	1,4 Kg
Air chamber diameter	57 mm
Piston travel	210 mm
Tubes	N° 2 lenght 1 m

QM22125A

QMX 125 ALUMINIUM  
DOUBLE ACT VACUUM  
PUMPT FOR VACUUM  
SYSTEMS



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>