

## 510 SELF

### Self loading trolley



510 Self is a self loading trolley to be used with stretcher Spencer 506 for rescue operation on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically.

#### Specific features

The coupling with the stretcher is simple and intuitive  
 Color coded commands  
 2 fixed wheels and 2 pivoting with large diameter to make movement easier  
 Brakes on rear wheels  
 Locking system for legs in closed position. This system allows to lift and move the trolley and stretcher assembly.

#### Technical data

Length	1980 ±10 mm
Width	590 ± 20 mm
Weight	34 ± 1 kg
Loading capacity	170 kg
Wheels	ø200 mm
Loading carriage wheel height	Min 49 ± 1cm Max 62 ± 1 cm
Frame material	Steel

Class I Med.Device according to Dir. 93/42/CEE

510 SELF TROLLEY SILVER ROLL-ON  
 ROLL-OFF FOR 506T,  
**ST10602C**



CND Classification	V0899
Registration number	159722

Rev.1 (04/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)