

ALTO T

Stretcher with hydraulic system



Spencer ALTO is a sturdy stretcher designed for transport of patients inside hospitals. It is characterized by an hydraulically operated height adjustment system that confers high versatility to the device and minimizes fatigue.

Specific Features

- Polyethylene board easy to clean
- Adjustable backrest
- Trendelenburg antishock position
- 4 pivoting wheels with brake 2 of which are semi-locked
- Foldable side bars *Safe bar*
- Equipped with basket for storage of tools or devices
- Oxygen Tank carrier installed on the main frame, suitable for 5/7/10 liters oxygen tanks
- Additional anterior oxygen tank carrier
- Foot operated hydraulic height adjustment system
- Can be shortened for a better handling in confined spaces
- Integrated foldable IV holder

Standard Equipment

- STX 702 two pieces black reflex metal belts

Technical data

Length	2000 ± 10mm
Length folded	1640 ± 20mm
Width	680 ± 10mm
Weight	72 ± 2 kg
Loading capacity	250 kg
Wheels	ø200 mm
Patient floor height	From 650 to 1060 ±20 mm
Basket dimensions	450x640x80 ± 10 mm
Tank container dim.	ø145x740 ± 10 mm (internal)
Anterior tank carrier dim.	ø110x340 ± 20mm (internal)
Frame material	Steel
Side bars length	680 mm
Side bars height	200 mm

Accessories

ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Class I Med.Device according to Dir. 93/42/CEE

ALTO T CHROME W/YELLOW
BOARD W/ HYDRAULIC SYSTEM

BO01041C

CND Classification V9099

Registration number 89557



Rev.0 (13/01/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it