

## Spencer 420

Transport chair with two wheels



Spencer 420 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

### Specific Characteristics

Easy to clean PVC sheets  
 2 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation  
 Telescopic front handles to facilitate the handling  
 Foldable to save space  
 Strong and light frame

### Standard equipment

ST00423A | Pair of belts for chair 420/425 orange

### Accessories

ST36451B | F 400 fixation system 20 g for chairs

### Technical data

Length opened w/closed handles	660 ± 10 mm
Length opened w/opened handles	1150 ± 10 mm
Width	510 ± 10 mm
Height opened	900 ± 10 mm
Thickness folded	210 ± 10 mm
Rear wheels	ø125 mm
Materials	Steel/Al/PU/PVC
Weight	10,5 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

STAIR CHAIR SILVER SPENCER 420S

ST10420A

CND Classification V9099

Registration number 93511

