

SPENCER 420

Transport chair with two wheels



The lightest chair is also characterised by the extreme simplicity of use and by being extremely compact in its closed position. This characteristic gives it extreme adaptability in terms of ambulance storage. The Spencer 420 chair is supplied with two polypropylene 50 mm [1.96 in] straps.

TECHNICAL DATA

High opened	900 mm
High folded	900 mm
Width	490 mm
Thickness folded	180 mm
Length	590 mm
Length with extended handles	1180mm
Weight	10,5Kg
Load capacity	150Kg
Seat and back material	Spentex®
Rear wheels	Ø125 x 32mm
Front wheels	Not present
Brake	Not present
Frame	Steel
Surface treatment	Painting
Footrest and telescopic handles	Alluminium

ST20420A

SPENCER 420
TRANSPORT CHAIR
WITH TWO WHEELS



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>