

Eco-B-Life

Adult disposable manual resuscitator



Disposable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The material of construction and the surface finish increase grip to the operator's hands
 Connection standard ISO 15/22 mm
 Low storage space required
 The kit includes rubber face mask size 4, overpressure valve and reservoir

Technical data

Materials	Rubber, PC
Bag volume	1600 ml
Stroke volume	770 ml with 1 hand 900 ml with 2 hands
Max. re expansion rhythm	45 bpm
Reservoir volume	2500 ml
Overpressure valve	60 cmH ₂ O

Optional accessories

BG02001A	Laxan transparent bag for manual resuscitator
----------	---

Class IIa med. dev. according to dir. 93/42/CEE

ECO B-LIFE - RUBBER ADULT RES.
 W/OVERP.V. MASK S.4

BG06009B

Registration number



57263

Rev.0 (09/05/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it