

## TANGO

Integrated paediatric and adult spine boards



In the world of emergency care, it is of great importance to optimize stocking space and to satisfy all requests for ergonomically advanced rescue equipment. Optimizing space and performance is the line Spencer drew on when introducing the new Tango. A new design that brings together solidity, dynamics, space and performance and at the same time sets a new standard for high-performance spine boards. Great performance and above all, great for its attention towards the patient; adult and child. An exclusive idea that allows two different profiles with one single product, reducing costs and stock volume. A paediatric and adult spine board always ready to collaborate together and each enhancing the other.

### TECHNICAL DATA

Dimensions:	445 x 55 x h1830 mm
Weight:	8 kg
Load capacity:	150 kg
Material:	LLPE
Number of handles:	14
X ray translucent :	Yes for X ray diagnosis
Radio compatible:	Yes, except for nuclear magnetic resonance
Floatation:	Yes (Weight of device)
Tapered angles:	Yes
Functioning temperature:	from -20 to +60 °C

Tango

ST02140C

