

T SERIES

T1



PVC bag suitable to store devices commonly used during intubation procedures.

Specific Characteristics

Semi transparent to allow view of the content
 Zip closure for easy access to the content
 2 female strap bands to allow attachment inside a bag equipped with male straps

Technical data

Dimensions	330x110x30 mm
Weight	50 g
Material	PVC

T1
 CB02070B

Rev.0 (10/02/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it