

PARAMEDIC

Complete medication bag



Professional bag equipped with main devices used for medication.

Specific features

Detachable semi-rigid internal structure
 Compartment with zip closure and transparent window between the join point of upper and lower part (approximate size of the compartment 280x60x80mm)
 3 pockets on the internal side of the upper flap, two of which are mesh type. Above one of the two mesh pockets, is sewn an elastic band to attach small accessories
 Removable fabric element with strap to fix small pockets
 Reflective inserts to increase visibility
 Handle and shoulder strap for easy transport

Equipment

Disposable sterile nitrile gloves (2pairs)
 TIV – automatic tourniquette
 Silk bandaid 2,5cm x 9,14m
 Strip bandages. 20pcs diff.sizes (2packs)
 Sterile compresses gauze 10x10cm 25pcs (4 packs)
 Burn S – Burn stop dressing 10x10cm
 Gauze bendage 5cm x 5m
 Gauze bendage 7cm x 5m
 Gauze bendage 10cm x 5m
 Spencer ICE PACK – instant ice
 Isothermal sheet 220x160 cm
 Kelly straight hemostatic forceps 143mm
 Scissors for bendages Lister 189mm
 Surgical forceps 140mm
 Jaw spreader

Flex mask – mouth to mouth mask with box
 Plastic pen light 140 mm
 Black phonendoscope
 Electronic thermometer
 Universal – Tuff cut scissors
 TWK033 – Aluminum torch
 DGX410 – Adult sphygmomanometer

Technical data

Dimensions	360x200x130 ± 10mm
Shoulder strap length	Max 125 cm
Empty bag weight	650 ± 20 g
Materials	Nylon, PVC

Med.dev. CE marked according to dir. 93/42/CE

PARAMEDIC M COMPLETE MEDICATION -
 EMERGENCY BAG
 CB07000C

Rev.0 (21/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it