

## PARAMEDIC

### Complete medication bag



Professional bag equipped with main devices used for medication.

#### Specific features

Detachable semi-rigid internal structure  
 Compartment with zip closure and transparent window between the join point of upper and lower part (approximate size of the compartment 280x60x80mm)  
 3 pockets on the internal side of the upper flap, two of which are mesh type. Above one of the two mesh pockets, is sewn an elastic band to attach small accessories  
 Removable fabric element with strap to fix small pockets  
 Reflective inserts to increase visibility  
 Handle and shoulder strap for easy transport

#### Equipment

Disposable sterile nitrile gloves (2pairs)  
 TIV – automatic tourniquette  
 Silk bandaid 2,5cm x 9,14m  
 Strip bandages. 20pcs diff.sizes (2packs)  
 Sterile compresses gauze 10x10cm 25pcs (4 packs)  
 Burn S – Burn stop dressing 10x10cm  
 Gauze bendage 5cm x 5m  
 Gauze bendage 7cm x 5m  
 Gauze bendage 10cm x 5m  
 Spencer ICE PACK – instant ice  
 Isothermal sheet 220x160 cm  
 Kelly straight hemostatic forceps 143mm  
 Scissors for bendages Lister 189mm  
 Surgical forceps 140mm  
 Jaw spreader

Flex mask – mouth to mouth mask with box  
 Plastic pen light 140 mm  
 Black phonendoscope  
 Electronic thermometer  
 Universal – Tuff cut scissors  
 TWK033 – Aluminum torch  
 DGX410 – Adult sphygmomanometer

#### Technical data

Dimensions	360x200x130 ± 10mm
Shoulder strap length	Max 125 cm
Empty bag weight	650 ± 20 g
Materials	Nylon, PVC

#### Med.dev. CE marked according to dir. 93/42/CE

PARAMEDIC M COMPLETE MEDICATION -  
 EMERGENCY BAG  
 CB07000C

Rev.0 (21/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)