

E-TUBE

Insulation tube for Energy



Insulation element to be connected to Energy heating bag to limit the temperature drop caused by thermal dissipation along the infusion tube.

Specific features

- Easy connection to Energy heating bag by means of straps
- Strap longitudinal closure to allow easy insertion of the infusion tube
- Padding to improve insulation of the infusion tube

Technical data

Dimensions	90x1500 ± 20mm
Weight	470 g
Materials	Nylon, PVC

Class I Med.Device according to Dir. 93/42/CEE

INSULATING TUBE FOR SPENCER
ENERGY

IF03031A

CND Classification



A0899

Registration number

304199

Rev.0 (23/01/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it