

# 10G HOLDER

For Inject S



Device to be installed inside ambulances to hold the infusion pump Inject S.

### Specific Characteristics

Extremely sturdy structure  
 High grip components to decrease vibrations and noise  
 Easy to fix by means of two holes  
 10G compliant according to EN 1789

### Technical data

Dimensions	185x110x205 ± 10mm
Weight	1,34 ± 0,1 kg
Materials	Steel

Class I Med.Device according to Dir. 93/42/CEE

SUPPORT FOR 10G PUMP  
 INJECT



Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)