

## Emergency

### Emergency bag



Practical bag with everything is needed for an efficient approach to primary resuscitation.

#### Specific Characteristics

High portability thanks to the shoulder strap  
 High visibility due to reflective bands  
 Compact and practical to use  
 Logical and ordered placing of the content thanks to the rigid separation of the main compartment  
 2 lateral pockets with zip closure  
 Reinforced bottom

#### Equipment

FG00605A	Magill adult forceps
DG02110B	Black phonendoscope
ST00630A	HT Isothermic sheet adult
IF00150C	IFX50 Tubular tourniquet
OX20130A	OXX30 O <sub>2</sub> anticrush tube
RM10892A	Cuffed mask mis.3
RM10893A	Cuffed mask mis.4
RM30700A	Moulder jaw spreader
OX10005E	O <sub>2</sub> bottle 0,5l
PR00304B	Pressure reducer
RM30400A	Young steel tongue holder
DG00603A	Tuff cut Universal scissors
DG01411C	Sphygmomanometer
RM10321A	Disposable non-sterile Guedel sz.1
RM10322A	Disposable non-sterile Guedel sz.2
RM10323A	Disposable non-sterile Guedel sz.3
RM10324A	Disposable non-sterile Guedel sz.4
RV71120B	RES-Q-VAC Manual suction device
BG02007A	Pro-B-Life adult PVC disposable resuscitator
DG08062D	DGX062 Pen Light blue no/batt

#### Technical data

Dimensions	490x230x220 mm
Bag material	Nylon
Bag weight	730 g
Kit weight	4,8 kg
Tank connection standard	UNI
Reducer outgoing pressure	3,5 bar
Reducer outgoing flow	6 l/min
Anticrush tube length	2,1 m
Suction device max vacuum	300 mmHg
Suction device Max flow	22 l/min
Manual resuscitator volume	2500 ml
Man. resusc. Overp. Valve calibration	60 mbar

EMERGENCY – COMPLETE  
 EMERGENCY BAG  
**CB08000A**



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)