


## TRACK 11

### Compact IV-bottle holder



Device intended to support and protect the containers of medicinal substances administered to the patient by venous infusion. Allows a significant reduction of the oscillations that may damage the containers of liquids.

#### Specific Characteristics

- Easy ceiling or recessed installation
- Robust and light structure
- 10g fixing of the bags
- Space saving device
- Can hold 250 or 500cc bags
- Practical to use thanks to the locking system and straps

#### Compliant with

EN 1789

#### Technical data

Dimensions closed	210x110x27 mm
Dimensions opened	50x120x180 mm
Materials	Acciaio, Al
Weight	290 g

Class I Med.Device according to Dir. 93/42/CEE

TRACK 11 - COMPACT IV-BOTTLE HOLDER IN ALUMINIUM

IF01061B

CND Classification V80

Registration number 805628



Rev.1 (01/12/2014)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)