

RESPONSE

Multi-use emergency bag orange



In the lower part of the soft bag, a series of elastic loops and a pocket that closes with Velcro®, allow for the placing of varied types of accessories.

On the fly-front opening of the upper part, there are elastic loops.

Technical details

Dimensions	465x270xh270 mm
Capacity	34 L
Weight	1,2 kg
Material	Oxford 840 D
Carry strap	Yes
N. of handles	1
Rigid double bottom	Yes
Bottom reinforcement	Yes
Reflective insert	Yes
N. of product fixation	21
Cylinder fixation	Yes
Int. material	Heat moulded plastic
Personalization	Yes, hot print 100x50 mm

Latex free

CB07010A

B.L.S. BAG RESPONSE
ORANGE