

RESPONSE

Multi-use emergency bag red



In the lower part of the soft bag, a series of elastic loops and a pocket that closes with Velcro®, allow for the placing of varied types of accessories.

On the fly-front opening of the upper part, there are elastic loops.

Technical details

| | |
|------------------------|--------------------------|
| Dimensions | 465x270xh270 mm |
| Capacity | 34 L |
| Weight | 1,2 kg |
| Material | Oxford 840 D |
| Carry strap | Yes |
| N. of handles | 1 |
| Rigid double bottom | Yes |
| Bottom reinforcement | Yes |
| Reflective insert | Yes |
| N. of product fixation | 21 |
| Cylinder fixation | Yes |
| Int. material | Heat moulded plastic |
| Personalization | Yes, hot print 100x50 mm |

Latex free

CB07011A

B.L.S. BAG RESPONSE
RED