

FIRST RESPONSE

Multi-use first response bag



Response is an ideal, multi-purpose bag for the paramedic or volunteer who needs a robust holder in which to organise all the equipment necessary for use in First-aid. Thanks to its compact structure and the variety of ways in which it can be used, it can be used on any occasion. It is made of tear proof Oxford nylon which is easily washable. It consists of a rigid under part and of a soft part. The lower part is in preformed shockproof plastic and is equipped with housing for oxygen tank. Both parts are easily accessible, thanks to the use of zip openings. A transparent, external pocket with access from inside, allows for the placing of lists or to personalize the bag in a rapid, safe and secure manner. Response is fitted with handles, adjustable strap which permits use as a bag, complete with spring catches and rings in acetate resin, which if necessary will permit belt attachment of the bag. In order to have the smaller aids conveniently at hand, the rigid compartment has, on the upper part, a large pocket with a zip, which offers protection to a series of elastic loops, thus reserving the lower part for the bulkier accessories.

TECHNICAL DATA

Dimensions	470 x 270 x h280 mm
Volume	35 lt
Weight	2,05 kg
Material	Nylon
Tank compartment	Si

CB07014A

FIRST RESPONSE MULTI-USE
FIRST RESPONSE BAG

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>