

BURN C

Fire blanket



This blanket can be used as protection for rescue workers entering or leaving a fire area. It can also be used to extinguish flames on a burn victim and, if necessary, be used as a large dressing, providing instant pain relief and beginning the healing process while victims are being transported to emergency care units.

Specific Characteristics

Sterile trauma blanket
Shoulder strap make it easy to carry
Suitable for: extinguish flame, large dressing, entry or exit from fire
No loose fibers
Relieves pain
Cools and provides the necessary water supply to the skin
Non-adherent and non-toxic

Technical data

Dimensions	160x120 cm
Idrogel composition	96% water, Melaleuca Alternifolia, emulsifiers
PH Gel	5,5-7
Sterilized by	Gamma ray
Blanket Material	Wool

Class IIb MD according to Dir. 93/42/CEE

BURN C 160 X 120 CM - FIRE BLANKET
BS04050



CND Classification	M040499
Registration number	67349

Distributed by Spencer Italia S.r.l.

Rev.0 (27/09/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it