

## MINI BURN KIT



Burn kit equipped with primary devices to be used for burns treatment.

### Specific features

- The devices are stored inside a semi-rigid case
- The handle allow easy transport
- Double zip closure to allow easy access to the content
- Two internal flaps with mesh pockets and straps to store accessories

### Equipment

- Spencer Ice Pack
- HT Adult isothermic sheet
- Sterile gauze 10x10 cm
- 2x BURN S 10 x 10 cm - burn stop dressing
- 10x BURN B - burn stop gel pack
- Gauze bandage 5m X 7cm
- Sterile disposable nitril gloves size L
- UNIVERSAL BLACK - tuff cut scissors

### Technical data

Case dimensions	280x230x85 ± 5mm
Overall weight	1,2 kg

Med.de. CE marked according to dir. 93/42/CEE

MINI BURN KIT - BURN KIT

**BS05005B**

Rev.0 (26/01/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)