

**MAX HP**

Case in ABS with internal compartments and wall support



The container is made in ABS, guarantee of a lightweight product with an exceptional resistance to shocks, high and low temperatures. It is composed of two specular compartments, complete with dividers to improve the placement and protection of the contents. To ensure individualization of these contents the two shells are protected and divided by transparent walls in Vitroflex. Neoprene trimming guarantees it is airtight. A simple system allows opening to 90 ° whilst still attached to its support. This support in ABS complete with wall bracket, allows the drawer to be safely positioned. The closing system is two levers which are activated with an intuitive rotation. The drawer is equipped with a convenient handle that makes transport immediate and secure.

**TECHNICAL DATA**

Dimensions:	443 x 147 x h338 mm
Weight:	8,5 kg
Material:	ABS
Wall fixation:	Yes
Closure:	2 Turning clips
Sterile gloves:	5
Face shield:	1
Physiological liquid sterile 500 ml:	3
Disinfectant bottle 500 ml:	2
Hydrophilic cotton box 100 g:	1
Tourniquet:	3
Silk plaster 2,5 cm:	2
Assorted plasters box 10 pcs:	2
Sterile gauze dressings 10 x 10 cm:	10
Sterile gauze dressings 18 x 40 cm:	2

Max HP

MX10100C