

MAX PRO

Drawer in polypropylene with internal compartments and wall bracket



The drawer, a single piece, is made in polypropylene. It is characterized by a top closure which presents important advantages during both the filling and use of this particular container.

The support is in plastic and the wall bracket ensures safe storage. The drawer is equipped with a convenient handle that renders transport immediate and secure.

TECHNICAL DATA

Dimensions:	250 x 90 x h190 mm
Weight:	1,4 kg
Material:	Polypropylene
Wall fixation:	Yes
Closure:	2 Turning clips
Sterile gloves:	2
Physiological liquid sterile 500 ml:	1
Hydrophilic cotton box 100 g:	1
Tourniquet:	1
Assorted plasters box 10 pcs:	1
Plaster roll 2,5 cm x 5 m:	1
Sterile gauze dressings 10 x 10 cm:	1
Sterile gauze dressings 18 x 40 cm:	1

Max Pro

MX10101C