

LIFE BAG 3

Emergency professional bag



Professional bag ideal to store rescue equipment or accessories, characterized by easy access thanks to the independent compartments.

Specific features

Equipped with handles and shoulder strap

Main compartment with adjustable separator and elastic bands to attach accessories

The lateral compartments have with zip closure and flaps that can be opened to 90° for an easy access to the content. Inside there are elastic bands to store vial or and accessories

Inside the lateral compartments, there are practical mesh pockets

The padding provides protection to the content

Technical data

Dimensions	400x260x210 mm
Volume	19 l
Weight	920 g
Materials	Nylon, PVC

LIFE BAG 3 - PROFESSIONAL EMERGENCY
BAG

CB07050A

Rev.0 (12/03/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it