

BLUE KIT

Kit sport



Compact bag equipped with devices frequently used by sport rescuers.

Included Equipment		Specific Characteristics	
BS04035D	BURN S (BA) 55 X 40 CM - BURN STOP DRESSING	Dual transport possibility thanks to the handles and shoulder strap	
CCF2000	12 SAFETY PINS	1 top pocket with zip closure	
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS	1 frontal pocket with zip closure	
DR10380K	GMX300 STERILE DISP. NYTRIL GLOVES Sz. L (1 PAIR)	2 side pockets with zip closure	
GM00300K	GMX300 -DISPOSABLE STERILE TWEEZERS	1 posterior pocket with strap closure	
GM00747K	GMX-747 TRIANGULAR BANDAGE TNT cm 96x96x136	Main compartment with zip closure and padding to protect the content	
GM02101K	GAUZE COMPRESSES STERILE 10X10CM 3PCS	Includes 5 colored bags with transparent upper side, fixed with strap	
GM02250K	CLOTH PLASTER 5CM X 5M	Internal divisor with 3 transparent pockets with zip closure and a big one with strap closure	
GM02405K	GAUZE BENDAGE 5M x5CM	1 inside pocket with zip closure in the internal side of the front flap	
GM02407K	GAUZE BANDAGE 5M X 7CM	Reflective bands to increase visibility	
GM02411K	GMX 401DISINFECTANT TOWEL 5PCS		
GM02510K	COTTON WOOL 100 gr.		
GM03011E	STERILE COMPRESSION BANDAGE 8X10 cm		
GM03215K	TAPE ADHESIVE PLASTER 1,25CM X 9M		
GM04021K	STRIP BANDAGES 20 PCS DIFFERENT SIZES 2PCS		
GM04100E	ELASTIC BRACE FOR MEDICAL 2pcs		
GM04242	AMUCHINA 500 ml.		
GM04250K	GMX 003 – DISINFECTANT 250 ml		
IF00426A	SPENCER ICE PACK PE BOX		
Technical data			
Dimensions		200x420x350 ± 10mm	
Internal bags dimensions		3 big 260x110x50 mm	
		2 small 130x100x40 mm	
Materiale		Nylon	
Volume		20 ± 0,5 L	
Weight of the empty bag		970 ± 20 g	
Medical devices CE marked according to dir. 93/42/CE			
BLUE KIT - SPORT KIT			
MX10305B			

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it