

UNIVERSAL

Tuff cut scissors



Tuff cuff scissors suitable to cut bandages, clothes and belts.

Specific features

Light and comfortable to hold

The terminal part is designed so as to be inserted between bandage and skin without causing injuries

Technical data

Dimensions	190x90x10 mm
Blade length	55 mm
Steam sterilizable at 134°	

Class I Med.Device according to Dir. 93/42/CEE

UNIVERSAL ORANGE - TUFF CUT SCISSORS

DG00601A

Rev.1 (12/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it