

**TOP 500**

Professional bag with integrated oxygen tank carrier



A very spacious and compact bag which can also accommodate an oxygen bottle.

- A large front pocket with zip and internal elastic loops
- Two mesh pockets with zip on internal flap
- Seven expanding mesh pockets on the back with elastic loops
- Two lateral pockets with zips
- A large central space which can accommodate a 2 or 5 L oxygen tank with elastic loops and elasticised mesh pockets
- Two pockets in elasticised mesh, with zips, inside the upper flap
- Removable and padded base, frontal and internal back parts
- Adjustable shoulder strap

**TECHNICAL DATA**

Dimensions	650 x 270 x h250 mm
Weight	1,5 kg
Capacity	42 L
Material	Nylon Ripstop

Top 500 CB09500B