


## TRI TAG



Triage cards

### Specific features

Ideal to deal with elevated numbers of victims requiring different treatment

Coloured on both front and backside, they have detachable, perforated and numbered strips with five different colours

### Technical data

Package dimensions <sup>(25sheets)</sup>	208x105x8 mm
Material	Paper

TRI TAG - TRIAGE CARDS

**TT02010B**

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)