


# Dolphin

## Universal linear life buoy



Linear life buoy to be used as buoyancy aid of victim during water rescue. Placed around the victim's abdomen, allows the floating in supine position.

### Specific features

- The segmented body confers easy grip and allow an easy placing around the victim's body
- The tow strap integrated in the main body, ensures high reliability
- The width of the shoulder strap tape, provide more comfort to the rescuer
- Metal carabiner and closure rings
- The handles offer great versatility
- The built-in face mask allows to carry out primary resuscitation procedures

### Standard equipment

TA09050A flex mask w/box blue logo spencer

### Optional

DY32401B HIGH VISIBILITY STORAGE BAG WITH REFLECTIVE TAPE

### Technical data

Main body dimensions	1200x175x85 ± 10 mm
Tow strap length	2200 mm
Belts width	30mm / 50mm shoulder strap
Materials	PE esp, Nylon
Weight	2,1 ± 0,1 kg
Minimum buoyancy force	100 N

Class I med. Dev. According to Dir. 93/42/CEE

DOLPHIN - SALVAGENTE  
LIN.UNIV.C/FLEX M.E MANIGLIE  
DY32500D



CND Classification V9099

Registration number 887374

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)