

TRY 744

Full Body manikin with electronics



Full body manikin studied for training in cardiopulmonary resuscitation (CPR) and external cardiac massage (MCE). The electronic part is equipped with led lights to indicate if the operations are properly carried out.

Specific Characteristics

Easy and fast assembly

Articulated body

Electronics shows:

Correct position of the hands

Correct cardiac compression

Correct ventilation

Electronics are compliant with EMC 89/336/EEC

Equipment

5 spare mouth-piece

5 spare airways with non rebreathing valve

4 Batteries type AA included

Technical data

Dimensions	127x53x30 cm
Weight	24 kg
Electronics compliant with EN60601-1-2	
Consolle specifications	
EN6100-4-2	6 kV Contact
EN6100-4-2	8 kV Air
EN61000-4-4	5kHz EFT/Burst

TRY 744FULL BODY-CPR UPPER

BODY MAN. WITH EL.P.

ST01250C

Distributed by Spencer Italia s.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it