

TRY 702

Realistic manikin



Full body manikin designed to provide effective training to victim handling, transportation, and extrication training. Used by Military forces, Fire, and Police Departments, Safety Teams, and Emergency Personnel for rescue and extrication.

Specific Characteristics

The fully articulated body, can be positioned in every conceivable situation and allow a realistic distribution of weights during recovery operations

Suitable for training in confined spaces, collapsed buildings, smoke rooms, and ladder carry-down protocols worldwide.

These manikins can be safely used in situations too hazardous or uncomfortable for human volunteers

Technical data

Height	1,65 m
Materials	PVC e galvanized steel
Weight	66 kg

TRY 702 - REALISTIC MANIKIN
ST01401C

Distributed by Spencer Italia s.r.l

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it