

## MATCH 3

### Multi-purpose scissors



1 Used to break glass.



2 Used to cut seatbelts.



3 Used to cut clothes.

Multi-purpose scissors with wall mount to provide safe installation inside any vehicle.

#### Specific features

- Interlocking fastening allows a quick access
- The support ensures that the device remains closed when not used
- The fastener is installed by means of 3 screws for maximum safety
- Rounded tip to limit risks of injury to the skin
- The shape of the plastic part covering the blade, allows to grasp the scissors using them as an evacuation hammer

#### Technical data

Scissors dimensions	247x105x18 mm
Dimensions with fastener	252x115x22 mm
Blade length	100 mm
Material	Steel, Nylon
Scissors weigh	374 g
Fastener weight	36 g

#### Spare parts

- GD02064 Wall fastener
- GD02069 Hooks for carpet fastening

Class I Med.Device according to Dir. 93/42/CEE

MATCH 3 MULTIFUNCTIONAL  
SCISSORS

**GD02060C**

**CND Classification** L01040199

**Registration number** 161376



NATO STOCK N° 6515-15-010-7023

Rev.3 (07/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)