

MATCH 3

Multi-purpose scissors



Used to break glass.



Used to cut seatbelts.



Used to cut clothes.

Multi-purpose scissors with wall mount to provide safe installation inside any vehicle.

Specific features

- Interlocking fastening allows a quick access
- The support ensures that the device remains closed when not used
- The fastener is installed by means of 3 screws for maximum safety
- Rounded tip to limit risks of injury to the skin
- The shape of the plastic part covering the blade, allows to grasp the scissors using them as an evacuation hammer

Spare parts

- GD02064 Wall fastener
- GD02069 Hooks for carpet fastening

Technical data

Scissors dimensions	247x105x18 mm
Dimensions with fastener	252x115x22 mm
Blade length	100 mm
Material	Steel, Nylon
Scissors weigh	374 g
Fastener weight	36 g

Class I Med.Device according to Dir. 93/42/CEE
MATCH 3 BLACK YELLOW HANDLE -
MULTI-PURPOSE SCISSORS

GD02065C

CND Classification	L01040199
Registration number	161406



Rev.3 (07/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it