

FORCE KIT 2

Force kit



Indicative image



Multifunctional kit that combines the main tools used to free rescue areas to make rescue operations quicker and safer.

Specific features

The compact bag with shoulder strap and internal rigid compartment, offer great portability and allows to store heavy tools.

Reflective bands increase visibility

Zip closure for quick access to the content

Frontal pocket with zip closure and a transparent pocket with access from the internal side of the bag

Technical data

Bag dimensions	480x265x135 mm
Weight	6,25 kg ± 0,2 kg

Content of the kit

Wind-proof torch with drops collector (2pcs)

Tow rope 5 m

Hatchet

Tongs

Crowbar

Leather gloves size 10, Cat. 1 (2 pairs)

Red warning light (2pcs)

Fireproof blanket EN1869 1,2x1,2 m

FORCE KIT 2 - FORCE KIT

GD04100B

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it