

## VISO BUMPER

Aneroid sphygmomanometer with large dial and protection



Wall installation sphygmomanometer, for the non –invasive measurement of arterial blood pressure.

### Specific features

- Dial with protection
- Spiral connection pipe to sphygmomanometer to avoid obstacles while using in confined spaces
- Dial with numbering at 20mmHg and division of 2 mmHg
- Compartment for housing the cuff and any small accessories
- Integrated elastic band in the compartment to prevent the leakage of the accessories
- Cuff with indications direction of application

### Conforms to

UNI EN ISO 81060-1

### Technical data

Extended cuff dimension	530x145 mm
Viso Bumper dimension	ø220x150 mm
Dial diameter	ø130 mm
Suitable for circumferences	from 25,4 to 40,6 cm
Pipe length	50 cm
Spiral pipe length	from 15 to 170 cm
Manometer range	0-300 mmHg
Material	PE, Nylon, rubber, Al
Total weight	1,1 ± 0,1 kg

MD class Im CE marked according to Dir. 93/42/CEE

VISO BUMPER – YELLOW  
ANEROID SPHYG. WITH LARGE  
DIAL AND PROTECTION



DG01600A

CND classification	C9006
Registration identifier	89767

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)