

VISO KOMPAK-T

Aneroid sphygmomanometer for table



Sphygmomanometer, for the non –invasive measurement of arterial blood pressure .

Specific features

- Large dial
- The wide base makes it very stable and the inclination allows easy reading
- Dial wit numbering at 20mmHg and division of 2 mmHg
- Cuff with indications direction of application

Technical data

Extended cuff dimension	530x145 mm
Sphygmom. Dimensions	166x176x138 mm
Suitable for circumferences	from 25,4 to 40,6 cm
Pipe length	50 cm
Manometer range	0-300 mmHg
Material	PE, Nylon, PC
Weight	0,8 ± 0,1 kg

Compliant with

UNI EN ISO 81060-1

MD class Im CE marked according to Dir. 93/42/CEE

VISO KOMPAK-T -LARGE FACE
ANEROID SPHYGMOM.FOR
TABLE



DG01502A

CND classification C9006

Registration number 89889

Rev.1 (25/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it