

## BLUE BAG 2

### Professional bag



Professional bag particularly suitable to store the main devices commonly used for measurement of vital signs.

#### Specific features

Frontal pocket with zip closure
Two side pockets with zip closure
Top pocket with zip closure
Main compartment with zip closure and equipped with:
- 5 internal bags of different colors with zip closure and transparent window
- Pocket with zip closure on the internal side of the frontal flap
- Internal separator characterized by a side made of tissue for the fixing of the internal bags; on the opposite side there are 3 transparent compartments with zip closure to store additional accessories
The main compartment is padded to provide protection to the content
Pocket on the back with strap closure
Maximum portability thanks to the handles and detachable shoulder strap
Reflective inserts to increase visibility

#### Technical data

Dimensions	390x210x340 ± 10mm
Volume	20 ± 0,5 l
Shoulder strap length	From 65 to 100cm
N° of bags	5
Bags color	Green, orange, red, yellow, blue
Materials	Nylon, PVC
Weight	1020 ± 20g

BLUE BAG 2 - PROFESSIONAL BAG  
CB09601C

Rev.0 (18/07/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)